

# SurCapt™

## Microbial Surface Detection Kit



The SurCapt™ Microbial Surface Detection Kit is an all-inclusive, ready-to-use kit which simplifies and quickens the process of environmental surface testing for biological contamination. The kit uses GreenLight® Technology for the rapid detection of microbial activity on cleanroom surfaces and the design reduces the likelihood of sample cross-contamination and mishandling of samples by the operator.

The SurCapt Kit is intended to be used with the GreenLight automated reader. The reader's internal LED sends light to the polymer sensor located in the kit vial. The optical system detects the decay of the fluorescence signal emission, referring to the oxygen depletion over time due to the microorganism's respiration. The system is flexible and allows either enumeration produced from an internal calibration, or a pass/fail test based on presence/absence of oxygen depletion over time. The reader software allows continuous measurements of oxygen depletion within the reader or endpoint measurements from externally incubated samples.

The SurCapt Kit with the Greenlight dedicated reader offers:

- A traditional, growth-based method with enhanced sensitivity due to the GreenLight Technology
- Non-destructive testing allows further analysis or speciation by traditional methods
- Suitability for Class A/ISO 5 classified areas and lower
- An ideal way to sample the most critical and hard-to-reach points identified by risk assessment
- High recovery of surface microbial contaminants using certified flocked swabs
- Faster time-to-results using an oxygen-sensitive probe assay
- Automated sample analysis eliminating human reading subjectivity
- Up to 24 samples can be automatically processed per carousel



## BENEFITS

### Time-to-Result

- Achieve results typically within 24 hours
- Increase sensitivity by maximizing recovery
- Advanced test tracking and reliability
- Easy-to-use automated reader with 21 CFR Part 11 compliance

### Cost Saving

- Easier disposal than agar plates
- No extra costs and risks due to reagent dispensing
- Reduces risks of false positives in aseptic production areas and related investigations

### Test Ready

- Fully prepared and disposable
- Minimizes operator handling
- Eliminates unintended contamination
- Each set of five samplers is triple-bagged and gamma irradiated
- H<sub>2</sub>O<sub>2</sub> - impermeable packaging

### Reduced Operator Handling

- Decreases time-consuming and operator-dependent manipulations
- Simplified standard operating procedures
- Faster sampling steps
- Easy lift-out interchangeable carousel

## FEATURES

- Standard TSB culture media
- FLOQSwab™ flocked swab with SRK™ for highest recovery swabbing (>70% each lot).
- Reduced variability on results
- Eliminate cleaning validation studies for culture media residues required for contact plate method
- Compact design and low maintenance reader
- Integrated incubator
- Intuitive screen monitor
- Easy-to-read dashboard with tabular graphics
- Result as Presence/Absence or as colony forming units (CFU)
- Secure results and 21 CFR part 11 compliance
- Compatible with Microsoft Windows®

## APPLICATIONS

- Microbiological surface sampling in cleanrooms and associated controlled environments
- Microbiological monitoring of equipment, operator clothing and gloves
- Easy monitoring of critical point in RABS and isolators

## Surface Detection Kit

<b>Materials</b>	Each individual gamma irradiated peel pouch kit contains: - SurCapt vial: polypropylene vial with integrated platinum porphyrin polymer and 10 ml of TSB - FLOQSwab: flocked swab in polypropylene tube with sponge and 1ml of Surface Rinse Kit™ (SRK)	
<b>Dimensions</b>	SurCapt vial - 15 ml capacity vial: - 120 mm length including red stopper - 17 mm diameter	Swab tube: - 175 mm total length including green stopper - 12 mm external diameter
<b>Filling Volume</b>	10 ml	
<b>Media</b>	Tryptic Soy Broth (TSB)	
<b>Additives</b>	Neutralizers included in the SRK buffer	
<b>Packaging</b>	Two bags containing five kits each are triple-bagged together, 20 of the triple-bagged units are placed inside a fourth plastic bag to produce a final package containing a total of 200 vial kits.	
<b>Storage Temperature</b>	5 °C - 25 °C (41 °F to 77 °F)	
<b>Approximate Detection</b>	10 <sup>0</sup> to 10 <sup>8</sup> cells	

## GreenLight Automated Reader

<b>Model</b>	Model 930-15 uses 15 ml capacity vials or adapters	<b>Operating Temperature</b>	18 °C to 24 °C (64.4 °F to 75.2 °F)
<b>Test capacity</b>	24 vials	<b>Storage Temperature</b>	10 °C to 35 °C (50 °F to 95 °F)
<b>Throughput</b>	Up to 216 tests in pass/fail mode	<b>Relative Humidity</b>	0 - 95% non-condensing
<b>Power (adaptor supplied)</b>	100-240 VAC, 50-60 Hz Output 72 W max	<b>Dimensions (H x W x D)</b>	31.8 x 30.50 x 38.14 cm (12.52 x 12.01 x 15.02 in)
<b>USB ports</b>	1 USB peripheral port 1 external Barcode Reader port	<b>Weight</b>	12.30 kg (27.12 lbs)
<b>Internal Incubation Temperature</b>	Ambient 2 °C to 40 °C (35.6 °F to 104 °F)	<b>Accessories</b>	Laptop computer (pre-loaded with GreenLight software), extra carousel (24 places), external Barcode Reader
<b>Language</b>	Manageable by end-user		

The SurCapt™ vial uses GreenLight® Technology, GreenLight® is a registered trademark of Luxcel BioSciences, Ltd.  
FLOQSwab™ and Surface Rinse Kit™ (SRK) are trademarks of COPAN Diagnostics Inc.  
Windows® is a registered trademark of Microsoft Inc.  
All other trademarks are the property of their respective company.  
Particle Measuring Systems, Inc. reserves the right to change specifications without notice.



### IRELAND

EMS Ireland

T: +353 (0)1 295 7373

E: [info@ems.ie](mailto:info@ems.ie)

W: [www.ems.ie](http://www.ems.ie)

### UNITED KINGDOM

EMS Particle Solutions

T: +44 01223 257704

E: [info@emsparticlesolutions.co.uk](mailto:info@emsparticlesolutions.co.uk)

W: [www.emsparticlesolutions.co.uk](http://www.emsparticlesolutions.co.uk)



### PMS HEADQUARTERS

5475 Airport Blvd

Boulder, Colorado 80301 USA

T: +1 303 443 7100, +1 800 238 1801

W: [www.pmeasuring.com](http://www.pmeasuring.com)