

SLS-20

Syringe Liquid Sampler

Specifications

Syringe size	1, 5, 10, 25, 50 (standard) ml
Sampling mode	Vacuum, Batch
Wetted surface materials	Kel-F®, PFA, PTFE, Kalrez®
Sample volume	Adjustable using software interface
Sample temperature	59 – 104 °F (15 – 40 °C)
Dimensions (l, w, h)	10.0 x 9.0 x 12.3 in (25.4 x 22.9 x 31.2 cm) not including particle counter
Weight	14 lb (6.4 kg) not including particle counter
Power	100 – 240 V, 50 – 60Hz
Environment	50 – 86 °F (10 – 30 °C) Humidity: Non-condensing



SOLUTIONS FOR YOUR INDUSTRY
Scan to learn more about our products!



HEADQUARTERS

5475 Airport Blvd
Boulder, Colorado 80301 USA
T: +1 303 443 7100, +1 800 238 1801

Instrument Service & Support
T: +1 800 557 6363

Customer Response Center
T: +1 877 475 3317
E: info@pmeasuring.com

www.pmeasuring.com
info@pmeasuring.com



EMS Particle Solutions Ltd
T: +44 (0) 1223 257 704
M: +44 (0) 7825 884292
E: info@emsparticlesolutions.co.uk

I EL

EMS (Ire) Ltd
T: +353 (0)1 295 7373
M: +353 (87) 254 5491
E: info@ems.ie

